

# Laboratory management software MetaBridge



## Discover the Brabender MetaBridge

The new software is characterized by its easy and intuitive handling. After log-in, the user finds all information about the device and a choice of options for his purpose on the start screen.

## The advantages

- Ready to use, no installation necessary
- User-friendly touch screen
- Test tracking independently from end device and location
- Responsive web design: automatically adjusted screen resolution
- Access through easy user log-in
- Touch support for tablets and smartphones
- Test tracking from multiple end devices at a time
- Security of tests and data: protection from unauthorized access
- Quick software updates available online

## Functions for measurement, evaluation and administration

Benefit from new and optimized functions:

- Administration mode for user access rights
- Central test administration
- Live tracking of tests with end time indication
- Record, evaluate, print and export test results
- Interactive editing of measuring data
- Automatic measuring data storage mode
- Integration of a reference curve
- Control, zoom and comment function for diagrams



The MetaBridge enables the operation of the Brabender instruments from a tablet or smartphone. Track your test results from your lab, your home office or on the go.



By a mouse click or touch, the user can choose one of the tiles that are easily distinguished by colour.



**Brabender® GmbH & Co. KG**

Kulturstr. 49-51 · 47055 Duisburg · Germany  
Phone: +49 203 7788-0  
metabridge@brabender.com  
www.brabender.com



Brabender agencies all over the world  
© 2019 Brabender® GmbH & Co. KG

All trademarks are registered.  
Subject to change of design  
and technical modification  
without notice.



## MetaBridge®

Smart Software for  
Interconnected Communication



... where quality is measured.

# Laboratory management software MetaBridge



The MetaBridge is an integrated software solution that interconnects your Brabender instruments with users and with each other. This enables users to access any instrument at any time and operate several instruments simultaneously as well.

As the software is web-based, it is platform-independent. This enables more than one user to log in at the same time and access the MetaBridge from anywhere around the world. A reading can be retrieved on a computer, tablet or smartphone from any location, provided the technical parameters are in place.

With the MetaBridge you are always interconnected.

The new Brabender Metabridge controller is ready for retrofitting your old devices. It enables you to upgrade your Brabender device for the Metabridge-Software.



More information:  
**Brabender MetaBridge**

## Functions

The web-based MetaBridge software solution enables you to monitor measurements being taken in the test lab on any instrument linked in to the MetaBridge live at any time.

- Easy
- Mobile
- Safe



## Data security

Access to the MetaBridge is granted via a password protected log-in. An effective rights management protects your measurement data.

- Authorization by password
- Automatic software updates
- Control by administrators

## Compatibility

One for all: all Brabender instruments can be interconnected with authorized users.

- For all with MetaBridge compatible Brabender devices
- Cross-linking between devices
- Consistency: same level of performance on all end devices (PC / Mac, smartphone, tablet)



## Operation

The MetaBridge was programmed focusing on intuitive user guidance and rapid handling learnability. Automatic default settings avoid errors.

- Individually adaptable
- Standard interface for all devices
- Automatic default settings



## Networking

Ready to communicate: MetaBridge enables you to monitor your readings live and record them simultaneously. The results can then be evaluated, printed and exported using the software.

- Always current data
- Central storage and management
- Simple feedback



## Ready-to-use anywhere

The MetaBridge does not require installation on a terminal device, but can be used immediately. Several users can access the MetaBridge simultaneously – even via the same log-in.

- Ready to use – fast and mobile
- Simultaneous access possible
- Remote Access by service staff

**... where quality is measured.**

MetaBridge®  
Brabender®